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Instructions for life science research use only. Not tested for use in diagnostic procedures. For *in vitro* use only.

Instructions For Use

VirSNIp SARS-CoV-2 Spike L452R

530

Cat.-No. 53-0793-96

Roche SAP n° 09 459 022 001

Kit with reagents for 96 PCR reactions 20 µl for genotyping of SARS-CoV-2 RNA [lyophilized]

1. Content, Storage and Expiry

Storage at Arrival:

1 Vial yellow cap 96 reactions CoV (lyophilized)

Store cooled or at ambient temperature
Do not freeze the lyophilized reagents.

- Kits are stable for one year after production (store 4°C to 25°C in the dark). See lot-specific expiry date.
• Reconstituted reagents are stable for two weeks if stored protected from light and cooled (2°C to 8°C).
• Dissolved reagent can be stored long-term if frozen (-15°C to -25°C). Avoid multiple freeze-thaw cycles.

2. Additional Reagents required

LightCycler® Multiplex RNA Virus Master
or 1-step RT polymerase

Cat.-No. 06 754 155 001
90-9999-96

3. Introduction

Hundreds of thousands of isolates have been sequenced since the SARS-CoV-2 genome was published (MN908947). The Spike L452R mutation is predicted to increase infectivity and host immune evasion and is a marker for the Californian CAL.20C variant; 81% GISAID entries are B.1.427 or B.1.429 lineage.

Table with 10 columns: Spike Prot. Variation, Genetic Variation, UK B.1.1.7, UK B.1.525, ZA B.1.351, Brazil P.1, Calif B.1.429, DK mink Clust V, Function, Effect, Single assay. Rows include mutations like del HV69/70, K417N, N439K, L452R, Y453F, E484K, N501Y, A570D, D614G, P681H, V1176F.

4. Description

A 111 bp long fragment is amplified and analyzed running a melting curve, using a 452R specific detection probe. The amplification of wild type 452L is not visible.

5. Specification

Sensitivity better than 500 copies viral RNA.

6. Sample Material and Extraction

Coronaviruses affect normally the lower respiratory system, but SARS-CoV-2 is found also in nose and throat. Typical clinical samples are throat and nasopharyngeal swabs, sputum, saliva or gargle solution. Product tested with heat-treated gargle solution. For RNA extraction see manufacturer's kit instructions.

7. Material Safety Data (MSDS)

This product is not hazardous (according to regulation (EC) No 1272/2008), not toxic, not IATA-restricted. Not from human, animal or plant origin. Product contains synthetic oligonucleotide primers and probes. According to OSHA 29CFR1910.1200, Australia [NOHSC:1005, 1008 (1999)] and EU Directives (EC) No 1907/2006 and (EC) No 2015/830 any products which do not contain more than 1% of a component classified as hazardous or classified as carcinogenic do not require a MSDS.

8. Instructions for Use

Instruction for Roche 480 instruments. Capillary LightCycler[®], LightCycler[®] 96, MyGo and BioRad CFX96 instruments give similar results. For other instruments use the SYBR Green melting option.

8.1. Programming Roche 480 Instruments (Standard ModularDx Program)

Detection Format 530 Channel **Set Quant Factor 10, Max Integration Time 1 sec**
 LightCycler[®] 480 Instrument: 483-533
 LightCycler[®] 480 II Instrument: 465-510
 cobas z 480 Analyzer (open channel): 465-510

Program Step:	RT Step	Denaturation	Cycling			Cooling
Parameter						
Analysis Mode	None	None	Quantification mode			None
Cycles	1	1	40-45			1
Target [°C]	55	95	95	60	72*	40
Hold [hh:mm:ss]	00:05:00	00:05:00	00:00:05	00:00:15	00:00:15	00:00:30
Ramp Rate [°C/s] 96	4.4	4.4	4.4	2.2	4.4	1.5
Ramp Rate [°C/s] 384	4.6	4.6	4.6	2.4	4.6	2.0
Acquisition Mode	None	None	None	Single	None	None

* 72°C step can be skipped. 95°C can be cut to 3 s, 60°C to 12 s. RT and Den to 3 min (total time 45 min) Table 1

8.1.1. Melting Analysis (may be added or programmed as second run)

Detection Format Hydrolysis Probe or SimpleProbe
 LightCycler[®] 480 Instrument: 483-533
 LightCycler[®] 480 II Instrument: 465-510
 cobas z 480 Analyzer (open channel): 465-510

Program Step:	Melting			Cooling
Parameter				
Analysis Mode	Melting Curves mode			None
Cycles	1			
Target [°C]	95	40	75	40
Hold [hh:mm:ss]	00:00:30	00:02:00	00:00:00	00:00:30
Ramp Rate [°C/s]	4.4	1,5	-	1.5
Acquisition Mode	-	-	Continuous	2.0
Acquisitions [per °C]	-	-	3**	None

Table 2

** Melting slope shall be 0.19 to 0.29°C per second. If reading more channels reduce the number of acquisitions/sec.

8.2. Experimental Protocol

- **Sample material:** Use aqueous nucleic acid preparations
- **Negative control:** Always run at least one no-template control (NTC) - replace the template NA with water.
- **Positive control:** Run a positive control

For an increased sensitivity use 10 µl nucleic acid per 20 µl reaction.
 Product tested for 10 µl reaction volume (192 reactions).

8.2.1. Preparation of Parameter-Specific Reagents (PSR, 96 reactions):

The reagent vial with the **yellow** cap contains the primers and probe to run 96+ PCR reactions.

Check for the orange pellet, then add 50 µl PCR-grade water, mix (vortex) and spin down.
 For robotic pipetting the volume can be extended to 55 µl (signals will decrease by 10-20%).

► Use 0.5 µl reagent per 20 µl PCR reaction.

8.2.2. Preparation of the Positive Control

- not provided

8.2.3. Preparation of the Reaction Mix

Multiply volumes by the number of reactions plus controls and one reserve, prepare in a cooled tube:

LightCycler® Multiplex RNA Virus Master		1-step RT Polymerase 90-9999-96	
10.4 µl	Water, PCR-grade	Water, PCR-grade	4.4 µl
0.5 µl	This reagent (PSR)	This reagent (PSR)	0.5 µl
4.0 µl	Roche Mastermix	TIB Mastermix	10.0 µl
0.1 µl	RT Enzym	-	-
15.0 µl		15.0 µl	

Table 3

Mix gently, spin down and transfer 15 µl (10 µl) per well.

Add 5 µl (10 µl) of sample or control to each well for a final reaction volume of 20 µl. Seal plate and centrifuge.

Start run

9. Typical Results


Figure 1. Green 452L (wild type), Red 453R, Purple 453F Left. Center: Melting Curve. Right 452L has a melting temperature of 50.3 (± 1)°C, 453R has a Tm of 56.8 (± 1)°C, 453F has a Tm of 45°C

10. Reading the Results

View results in the 530 channel. Use 'Tm calling'. The negative control (NTC) must show no signal.

Channel 530 Amplification	Channel 530 Melting analysis	Channel 530 NTC Control	Result
Not relevant	Not relevant	Negative / no peak	No virus amplified / not detectable
Invisible	Tm ~ 45°C*	Negative	SARS Spike 453F
Invisible	Tm ~ 50°C*	Negative	SARS Spike 452L
Invisible or low	Tm ~ 57°C*	Negative	SARS Spike 452R
Not relevant	Not relevant	Positive	Contamination Repeat experiment

Tm values shift depending on the instrument, speed of heating, mastermix, salt contents and detection format.

* Temperatures with 1step RT pol. 90-9999-96 are 3-4°C higher.

Single peak with lower Tm values are an indication for the presence of another mutation in the probe region.

The assay detects the genetic situation and not a strain; the correlation to a reference strain describes the most likely assignment for European isolates isolated winter 2020/2021.

11. References

Genomic characterisation of emergent SARS-CoV-2 lineage in UK defined by novel set of spike mutations. Rambaut et al., 2020 www.ecdc.europa.eu/en/publications-data/threat-assessment-brief-rapid-increase-sars-cov-2-variant-united-kingdom
 Mutations in SARS-CoV-2 spike protein and RNA pol. are associated with COVID-19 mortality risk. Hahn et al., 2020 www.rki.de/DE/Content/InfAZ/N/Neuartiges_Coronavirus/DESH/Bericht_VOC_05022021.pdf?__blob=publicationFile

12. Multiplex PCR Compatibility

This SNP assay can be combined with 51-0776-96 SARS E+N and either EAV spiked extraction control or UBC human mRNA extraction control or with the complete kit 60-0770-96 Sarbecovirus E+N+UBC.

Multiplex PCR and Instrument Compatibility

Color Comp 40-0320 mandatory only for Multiplex PCR using more channels

500	530	580	610	640	660
	SNP				
SarbecoV	SNP	UBC mRNA			
SarbecoV	SNP				UBC
SarbecoV	SNP				EAV
	SNP	SARS N	SARS E		UBC

480 II	z 480	LC96	LC2.0	Nano
X	X	X	X	X
X	X	X		
X	X	X		
X	X	X		
X	X	X		

Table 3

13. Version History

V210313 Release version

2020-03-13

Certificate of Analysis (CoA)					
Lot n° 5025 Expiry : YYYY-MM-DD					
Tm range	L452 48-51°C	452R 56-58°C	453F 44-46°C	Cp range	PC -
Measured					passed ✓
Signal level	2-10	2-10	2-10		✓
Measured					✓
Negatives	10/10				✓
Note: Cp (crossing point) values collected with pDNA (single target PCR). Fluorescence (FL) levels depend on instrument settings and may vary. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relative constant (ΔCp).					
DOM (manufactured): YYYY-MM-DD			QC Acceptance: YYYY-MM-DD		
We, the undersigned, certify that the product designated above has been obtained in accordance with the rules of production and quality control.					
Name(s) :					
<i>Name1</i>			<i>Name2</i>		

TIB MOLBIOL Syntheselabor GmbH | Eresburgstr. 22-23 | D-12103 Berlin | Germany
 Tel. +49 30 78 79 94 55 | FAX +49 78 79 94 99 | dna@tib-molbiol.de | WWW.TIB-MOLBIOL.COM
 Geschäftsführer (CEO): Olfert Landt | Register HRB 93163 B | Registergericht Berlin Charlottenburg

Instructions For Use

VirSNIp SARS-CoV-2 Spike P681R

530

Cat.-No. 53-0806-96

Roche SAP n° 09 496 742 001

Kit with reagents for 96 PCR reactions 20 µl for genotyping of SARS-CoV-2 RNA [lyophilized]

1. Content, Storage and Expiry

Storage at Arrival:

1 Vial yellow cap 96 reactions SARS CoV (lyophilized)

Store cooled or at ambient temperature Do not freeze the lyophilized reagents.

- Kits are stable for one year after production (store 4°C to 25°C in the dark). See lot-specific expiry date.
• Reconstituted reagents are stable for two weeks if stored protected from light and cooled (2°C to 8°C).
• Dissolved reagent can be stored long-term if frozen (-15°C to -25°C). Avoid multiple freeze-thaw cycles.

2. Additional Reagents required

LightCycler® Multiplex RNA Virus Master or 1-step RT polymerase

Cat.-No. 06 754 155 001 90-9999-96

3. Introduction

The SARS-CoV-2 genome was published 11.1.20 (Genbank MN908947). Hundreds of thousand isolates have been sequenced. The Indian origin variant B.1.617 carries the mutations 452R, 484Q and 681R :

Table with 10 columns: Spike Prot. Variation, Genetic Variation, UK B.1.1.7, UK B.1.525, ZA B.1.351, Brazil P.1, Calif B.1.429, India B.1.617, Function, Effect, Single assay. Rows include mutations like del HV69/70, K417N, N439K, L452R, Y453F, E484K/Q, N501Y, A570D, D614G, P681H/R, V1176F.

4. Description

A 103 bp PCR long fragment is amplified and analyzed with a melting curve using a 681R-specific probe. The amplification of 681H and 681P is not visible.

5. Specification

Sensitivity better than 100 copies viral RNA.

6. Sample Material and Extraction

Coronaviruses affect normally the lower respiratory system, but SARS-CoV-2 is found also in nose and throat. Typical clinical samples are throat and nasopharyngeal swabs, sputum, saliva or gargle solution. Product tested with heat-treated gargle solution. For RNA extraction see manufacturer's kit instructions.

7. Material Safety Data (MSDS)

This product is not hazardous (according to regulation (EC) No 1272/2008), not toxic, not IATA-restricted. Not from human, animal or plant origin. Product contains synthetic oligonucleotide primers and probes.

According to OSHA 29CFR1910.1200, Australia [NOHSC:1005, 1008 (1999)] and EU Directives (EC) No 1907/2006 and (EC) No 2015/830 any products which do not contain more than 1% of a component classified as hazardous or classified as carcinogenic do not require a MSDS.

8. Instructions for Use

Instruction for Roche 480 instruments. Capillary LightCycler[®], LightCycler[®] 96, MyGo and BioRad CFX96 instruments give similar results (FAM channel). For other instruments use SYBR Green melting option.

8.1. Programming Roche 480 Instruments (Standard ModularDx Program)

Detection Format 530 Channel **Set Quant Factor 10, Max Integration Time 1 sec**
 LightCycler[®] 480 Instrument: 483-533
 LightCycler[®] 480 II Instrument: 465-510
 cobas z 480 Analyzer (open channel): 465-510

Program Step:	RT Step	Denaturation	Cycling			Cooling
Parameter						
Analysis Mode	None	None	Quantification mode			None
Cycles	1	1	40-45			1
Target [°C]	55	95	95	60	72*	40
Hold [hh:mm:ss]	00:05:00	00:05:00	00:00:05	00:00:15	00:00:15	00:00:30
Ramp Rate [°C/s] 96	4.4	4.4	4.4	2.2	4.4	1.5
Ramp Rate [°C/s] 384	4.6	4.6	4.6	2.4	4.6	2.0
Acquisition Mode	None	None	None	Single	None	None

* 72°C step can be skipped. 95°C can be cut to 3 s, 60°C to 12 s. RT and Den to 3 min (total time 45 min) Table 1

8.1.1. Melting Analysis (may be added or programmed as second run)

Detection Format Hydrolysis Probe or SimpleProbe
 LightCycler[®] 480 Instrument: 483-533
 LightCycler[®] 480 II Instrument: 465-510
 cobas z 480 Analyzer (open channel): 465-510

Program Step:	Melting			Cooling
Parameter				
Analysis Mode	Melting Curves mode			None
Cycles	1			
Target [°C]	95	40	75	40
Hold [hh:mm:ss]	00:00:30	00:02:00	00:00:00	00:00:30
Ramp Rate [°C/s]	4.4	1,5	-	1.5
Acquisition Mode	-	-	Continuous	
Acquisitions [per °C]	-	-	3**	None

Table 2

** Melting slope shall be 0.19 to 0.29°C per second. If reading more channels reduce the number of acquisitions/sec.

8.2. Experimental Protocol

- **Sample material:** Use aqueous nucleic acid preparations
- **Negative control:** Always run at least one no-template control (NTC) - replace the template NA with water.
- **Positive control:** Run a positive control - replace the template NA with the provided Positive Control.

For an increased sensitivity use 10 µl nucleic acid per 20 µl reaction.
 Product tested for 10 µl reaction volume (192 reactions).

8.2.1. Preparation of Parameter-Specific Reagents (PSR, 96 reactions):

The reagent vial with the **yellow** cap contains the primers and probe to run 96+ PCR reactions.

Check for the orange pellet, then add 50 µl PCR-grade water, mix (vortex) and spin down.
 For robotic pipetting the volume can be extended to 55 µl (signals will decrease by 10-20%).

► Use 0.5 µl reagent per 20 µl PCR reaction.

8.2.2. Preparation of the Positive Control

- not provided -

8.2.3. Preparation of the Reaction Mix

Multiply volumes by the number of reactions plus controls and one reserve, prepare in a cooled tube:

LightCycler® Multiplex RNA Virus Master		1-step RT Polymerase 90-9999-96	
10.4 µl	Water, PCR-grade	Water, PCR-grade	4.4 µl
0.5 µl	This reagent (PSR)	This reagent (PSR)	0.5 µl
4.0 µl	Roche Mastermix	TIB Mastermix	10.0 µl
0.1 µl	RT Enzym	-	-
15.0 µl		15.0 µl	

Table 3

Mix gently, spin down and transfer 15 µl (10 µl) per well.

Add 5 µl (10 µl) of sample or control to each well for a final reaction volume of 20 µl. Seal plate and centrifuge.

Start run

9. Typical Results


Figure 1. Green 681H, blue 681P, red 681R, grey NTC Left Only 681R visible Center 681H and 681P can be not distinguished, 681R has a T_m of $64.5 (\pm 2)^\circ\text{C}$. Right Green 484E 681P (wild type/ B.1.1.7/ B.1.351 lila P.1), red 484Q 681R (India), yellow 484Q 681P (B.1.177), light blue 484E 681 other mutation (untypical B.1.351).

10. Reading the Results

Use the Second Derivative Maximum method (Automated (F" max). View results in the 530 channel. The negative control (NTC) must show no signal. For the melting curve analysis use 'Tm calling'.

Channel 530 Amplification	Channel 530 Melting analysis	Channel 530 NTC Control	Result
Not relevant	Not relevant	Negative / no peak	No virus amplified / not detectable
Invisible	$T_m \sim 55.0^\circ\text{C}^*$	Negative	SARS Spike 681H
Invisible	$T_m \sim 55.5^\circ\text{C}^*$	Negative	SARS Spike 681P
Visible	$T_m \sim 64.5^\circ\text{C}^*$	Negative	SARS Spike 484Q (eg. B.1.617)
Not relevant	Not relevant	Positive	Contamination Repeat experiment

T_m values shift depending on the instrument, speed of heating, mastermix, salt contents and detection format.

* Temperatures with 1step RT pol. 90-9999-96 are 3-4°C higher.

Single peak with lower T_m values are an indication for the presence of another mutation in the probe region.

The assay detects the genetic situation and not a strain; the correlation to a reference strain describes the most likely assignment for European isolates isolated winter 2020/2021.

11. References

Genomic characterisation of emergent SARS-CoV-2 lineage in UK defined by novel set of spike mutations. Rambaut et al., 2020 www.ecdc.europa.eu/en/publications-data/threat-assessment-brief-rapid-increase-sars-cov-2-variant-united-kingdom
 Circulating SARS-CoV-2 spike var. N439K maintains fitness while evading antibody-mediated immunity. Thomson et al., 2020
 Mutations in SARS-CoV-2 spike protein and RNA pol. are associated with COVID-19 mortality risk. Hahn et al., 2020
www.rki.de/DE/Content/InfAZ/N/Neuartiges_Coronavirus/DESH/Bericht_VOC_05022021.pdf?__blob=publicationFile

12. Multiplex PCR Compatibility

This SNP assay can be combined with 51-0776-96 SARS E+N and either EAV spiked extraction control or UBC human mRNA extraction control or with the complete kit 60-0770-96 Sarbecovirus E+N+UBC.

Multiplex PCR and Instrument Compatibility

Color Comp 40-0320 mandatory only for Multiplex PCR using more channels

500	530	580	610	640	660
	SNP				
SarbecoV	SNP	UBC mRNA			
SarbecoV	SNP				UBC
SarbecoV	SNP				EAV
	SNP	SARS N	SARS E		UBC

480 II	z 480	LC96	LC2.0	Nano
X	X	X	X	X
X	X	X		
X	X	X		
X	X	X		
X	X	X		

Table 3

53-0806 VirSNIp 681R and 53-0808 VirSNIp 484Q can be mixed.

13. Version History

V210404 Release version

2021-04-28

Certificate of Analysis (CoA)						
Lot n° 5075 Expiry : YYYY-MM-DD						
Tm range Measured	681H	681P	681R	Cp range	PC	passed
	54-57°C	54-57°C	64-66°C		-	✓
Signal level Measured	2-10	2-10	2-10		-	✓
Negatives	10/10					✓
Note: Cp (crossing point) values collected with pDNA (single target PCR). Fluorescence (FL) levels depend on instrument settings and may vary. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relative constant (Δ Cp).						
DOM (manufactured): YYYY-MM-DD			QC Acceptance: YYYY-MM-DD			
We, the undersigned, certify that the product designated above has been obtained in accordance with the rules of production and quality control.						
Name(s) :						
<i>Name1</i>			<i>Name2</i>			

TIB MOLBIOL Syntheselabor GmbH | Eresburgstr. 22-23 | D-12103 Berlin | Germany
 Tel. +49 30 78 79 94 55 | FAX +49 78 79 94 99 | dna@tib-molbiol.de | WWW.TIB-MOLBIOL.COM
 Geschäftsführer (CEO): Olfert Landt | Register HRB 93163 B | Registergericht Berlin Charlottenburg

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